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※ Product code : 11321 (20Tests/Kit)

INTENDED PURPOSE

BioTracer™ Influenza A&B Test is intended for the qualitative determination of influenza A virus and influenza B virus from nasal / throat swab specimens of symptomatic patients.

EXPLANATION OF THE TEST

Human influenza A and B viruses are major infectious pathogens of the respiratory tract that affect people of all ages, with the highest rates of mortality and severe morbidity in the elderly and very young. Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hrs of onset of illness. It is important to distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction(PCR), serology, and viral culture.

BioTracer™ Influenza A&B Test is a qualitative rapid test to detect influenza A or B antigens from respiratory specimens of symptomatic patients

PRINCIPLE OF THE METHOD

BioTracer™ Influenza A&B Test is based on the principle of an immunochromatography in vitro test for the qualitative determination of influenza A or influenza B. When the test strip is put into the test tube including specimens dissolved in extraction solution, the specimens moves to the conjugate pad via sample pad and resuspends anti-influenza A or anti-influenza B conjugated gold complex that is dried on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with anti-influenza A or anti-influenza B that is immobilized on membrane. If influenza A or influenza B is present enough in the sample, a colored band is appeared. If there is no influenza A or influenza B in the sample or not sufficient, the area will remain colorless. As the control band(C) which is red or purple color is appeared, it means the test is working properly and the result is valid.

CONTENTS

1. Test strip
2. Specimen extraction tube
3. Specimen collection swab
4. Influenza A positive control swab
5. Influenza B positive control swab
6. Influenza Negative control swab
7. Package insert

SPECIMEN COLLECTION AND PREPARATION

1. Specimen collection

Nasal swab

Throat swab

- 1) Nasal swab specimens
To collect a nasal swab specimens, insert the sterile swab into the nostril that presents the most secretion under visual inspection. Using gentle rotation push the swab until resistance is met at the level of the turbinate(less than one inch into the nostril). Rotate the swab a few times against the nasal wall.
- 2) Throat swab specimens
To collect a throat swab specimens, rub the sterile swab on both tonsillar surfaces and the posterior pharynx

2. Specimen transport and storage

Testing should be performed immediately after the specimens have been collected. The specimens extracted from swab may be stored at 1~8°C for up to one day. For long term storage, specimens should be kept below -20°C.

TEST PROCEDURE

Allow the specimens and sealed pouch containing the test strip to room temperature prior to use.

▶ Specimen extraction


All swab samples

- 1) Insert the patient swab sample into the test tube containing extraction solution after removing the cap. Swirl the swab at least 3 times while pressing the head against the bottom and side of test tube.
- 2) Swirl the swab head against the inside of the test tube after rolling the swab head at least 5 times when you remove it.

▶ Reaction with test strip

1. Remove the test strip from the foil pouch.
2. Place the test strip into the test tube in the direction of the arrows on the test strip pointing down.
3. Interpret test result at 15 minutes. Do not interpret the test result after 15 minutes.

※Caution : Perform the test immediately after removing the test strip from the foil pouch.


READING AND INTERPRETATION OF RESULT

Use the 'Reference paper for interpretation' provided to confirm the position of test band('T1', 'T2' band) and control band('C' band) of the test strip after performing test.

1. Control band(C) means that the test is working properly.
2. Test band(T1, T2) indicates the test results.

NEGATIVE :

The presence of only control band indicates a negative result.

POSITIVE :

- 1) Positive for Influenza virus A:

The presence of two color bands ("T1" band and "C" band) indicates a positive for Influenza A.


- 2) Positive for Influenza virus B :

The presence of two color bands ("T2" band and "C" band) indicates a positive for Influenza B.


- 3) Positive for Influenza virus A/B:

The presence of three color bands ("T1", "T2" band and "C" band) indicates a positive for Influenza A/B.

INVALID :

If control band('C' band) is not appeared after performing test, the result is considered invalid. The directions may not have been followed correctly or the strip may have been deteriorated. It is recommended that the specimen be re-tested with the new strip.

FOLLOW-UP ACTION

1. **BioTracer™ Influenza A&B Test** is a qualitative test kit. For knowing the quantitative concentration of Influenza, it should be performed with additional quantitative test.
2. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

LIMITATIONS OF THE METHOD

1. **BioTracer™ Influenza A&B Test** is for the qualitative detection of influenza A and B viral antigens. The test performance depends on antigen load and may not correlate with cell culture performed on the same specimen. Negative test results are not intended to rule out other non-influenza viral infections.
2. Sensitivity can differ with various strains of influenza due to difference in antigen expression. Specimens might contain new, non-identified strains of influenza that express varying amounts of antigen.
3. As with all diagnostic assays, the results obtained with this test kit yield data that must be used only as an adjunct to other information available to the physician.
4. Individuals who received nasally administered influenza vaccine may have positive test results for up to three days after vaccination.
5. Monoclonal antibodies may fail to detect(or with less sensitivity) influenza A viruses that have undergone minor amino acid changes in the target epitope region.

PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

A study was performed using 400 positive and negative specimens. Each specimen was assayed with the **BioTracer™ Influenza A&B Test** and a commercially available DFA after cell culture according to the respective instructions for use. Compared to culture and confirmed for influenza A or influenza B by DFA, the **BioTracer™ Influenza A&B Test** correctly identified 95.5%(191/200) positive specimens and 99.9%(200/200) negative specimens.

		BioTracer™ Influenza A&B Rapid Test			Total
		A	B	Negative	
Cell Culture Test	A	97	—	5	102
	B	—	94	4	98
	Neg.	—	—	200	200
Total		97	94	209	400

※ Sensitivity : Influenza A 95.10%(97/102), Influenza B 95.92%(94/98)

Specificity : more than 99.9%(200/200)

2. Analytical specificity and cross reactivity : **BioTracer™ Influenza A&B Test** was evaluated with 8 bacterial and 8 viral isolates. These studies had been performed to demonstrate that have been no cross reactivity with 8 bacterial panels and 8 virus panels on the **BioTracer™ Influenza A&B Test**.

3. Analytical sensitivity : Analytical sensitivity was established using a total 12 strains of influenza virus : 7 influenza A and 5 influenza B

4. Precision

- 1). Within run performance test was determined by one analyst with ten strips each of 3 different lots containing different concentration every influenza type. There was no variation between the test strips of each lot and between three different lots.
- 2) Between run performance test was determined by three analysts with 3 different lots containing different concentration every influenza type. There was no variation between different analysts.

5. Analytical Interference test

The following potential interfering substances were tested and were found to have no effect on the performance of the **BioTracer™ Influenza A&B Test**.

Potential interfering substances	Concentration
Acetyl Salicylic Acid	20 mg/mL
Acetamidophenol	10 mg/mL
Chlorpheniramine	5 mg/mL
Diphenhydramine HCl	5 mg/mL
Guaiacol Glyceryl Ether	20 mg/mL
Oxymetazoline HCl	10 mg/mL
Phenylephrine HCl	100 mg/mL
Fexofenadine	500mg/mL
Amantadine	500mg/mL
Ribavirin	500mg/mL
Pseudoephedrine HCl	20mg/mL
Ibuprofen	10mg/mL

WARNING AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Do not touch the membrane of the test strip by hands or it may affect the test result.
3. Treat all specimens as being potentially infectious. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
4. Any waste used for testing should be autoclaved at 121°C for 15minutes before disposal or dispose according to laboratory regulations for infectious waste.
5. This product is intended for qualitative test. For precise concentration of the specimen, quantitative assay is needed.
6. This test kit contains a little of sodium azide. Avoid contact with skin, eyes.
If you get the solution on your body, immediately wash thoroughly with flushing water. And contact a doctor if you need to.
7. Do not use test kit beyond the expiration date.
8. Avoid moisture, direct sunlight and heat, and keep it in room temperature.
9. The test device is sensitive to humidity. Pull out the test strip from the foil pouch right before use.
10. This test kit should be used as the first line diagnosis only.
Since the test may give you the false positive or false negative result, final diagnosis should not be made only by this product's test result.
11. The confirmation test would be necessary. And the final diagnosis should be made by doctor and other clinical data in conjunction with inquiries.
12. Do not eat or smoke while handling specimens.
13. Avoid splashing or aerosol formation.
14. Do not re-use the test strip.

※ Use the strip with the extraction tube in same box.

* Cautions for the specimen collection swab

- (1) Do not use the swab if it is visibly damaged or its package is open.
- (2) Do not re-use and re-sterilize the swab.
- (3) Do not use excessive force when collecting samples from patients as this may result in accidental breakage of the swab shaft.
- (4) Do not use it beyond the expiration date.
- (5) Keep it in room temperature where also dry and dark.
- (6) Use it only for the intended purpose.
- (7) The swab used for testing may contain infectious micro-organism; After use, the swab must be disposed as infectious waste.

STORAGE AND SHELF-LIFE

BioTracer™ Influenza A&B Test should be stored at 2-30 °C (36-86 °F). The test strip is sensitive to humidity as well as to heat. Do not use it beyond the expiration date, 24 months from manufacturing.

LITERATURE REFERENCES

1. Barenfanger, J., C. Drake, T. Mueller, T. Troutt, J. O'Brien, and K. Guttman. R-Mix cells are faster, at least as sensitive and marginally more costly than conventional cell lines for the detection of respiratory viruses. J. Clin. Virol. 2001; 22:101-110.
2. Glezen W. Emerging infections: pandemic influenza. Epidemiol Rev 1996;18:64-72.

Manufactured by
NanoEnTek Inc.

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,
Gyeonggi-do, 18531, Korea
Tel : +82-2-6220-7942 / Fax : +82-2-6220-7999
ivdst@nanoentek.com
www.nanoentek.com

EC REP

MT Promedt Consulting GmbH
Altenhofstrasse 80, 66386
St. Ingbert, Germany